



Sampling Line

Features

The sampling line are intended to interface with the patient via nasal cannula, nares connector or standard sampling port connectors

Environment of use- hospital and sub-acute settings.

Patient population-non-intubated patients who are breathing spontaneously.

Introduction for Use

- Remove the sampling line from the outer wrapper and connect to the patient breathing circuit
- Connect the sampling line to the monitor
- Do not use the sampling line for venous, arterial or enteral access

Specifications

Sampling Line Tube: ID1.0mm * OD2.2mm

Caution

Do not re-use ,do not use if packaged is damaged

Warnings

Do not attempt to clean disinfect ,or sterilizer the devices as it could lead contamination of instrument and may cause infection to patients and users